

## Introducer Sheath [INSTRUCTION FOR USE]

### 1. INTENDED PURPOSE

The product is intended to enlarge the incision of percutaneous intervention, which assists catheter into blood vessel by professional physician.

### 2. CONTRAINDICATION

No absolute contraindications when used correctly by medical personnel.

### 3. HANDING AND STORAGE

Stored in the relative humidity of not more than 80%, non-corrosive gases, a common, dry, well-ventilated and clean environment.

### 4. STRUCTURE AND SPECIFICATIONS

The product is consist of introducer sheath and dilator.

Specifications:

Specifications	Introducer sheath Size (Fr.) x Length (cm)	Dilator Size (Fr.) x Length (cm)	Specifications	Introducer sheath Size (Fr.) x Length (cm)	Dilator Size (Fr.) x Length (cm)
BR-RA0407	4F X 7cm	4F X 13cm	BR-RA0723	7F X 23cm	7F X 29cm
BR-RA0507	5F X 7cm	5F X 13cm	BR-FA0409	4F X 9cm	4F X 15cm
BR-RA0607	6F X 7cm	6F X 13cm	BR-FA0511	5F X 11cm	5F X 18cm
BR-RA0707	7F X 7cm	7F X 13cm	BR-FA0611	6F X 11cm	6F X 18cm
BR-RA0411	4F X 11cm	4F X 18cm	BR-FA0711	7F X 11cm	7F X 18cm
BR-RA0511	5F X 11cm	5F X 18cm	BR-FA0811	8F X 11cm	8F X 18cm
BR-RA0611	6F X 11cm	6F X 18cm	BR-FA0911	9F X 11cm	9F X 18cm
BR-RA0711	7F X 11cm	7F X 18cm	BR-FA1211	12F X 11cm	12F X 18cm
BR-RA0416	4F X 16cm	4F X 23cm	BR-FA0419	4F X 19cm	4F X 26cm
BR-RA0516	5F X 16cm	5F X 23cm	BR-FA0523	5F X 23cm	5F X 29cm
BR-RA0616	6F X 16cm	6F X 23cm	BR-FA0623	6F X 23cm	6F X 29cm
BR-RA0716	7F X 16cm	7F X 23cm	BR-FA0723	7F X 23cm	7F X 29cm
BR-RA0423	4F X 23cm	4F X 29cm	BR-FA0823	8F X 23cm	8F X 29cm
BR-RA0523	5F X 23cm	5F X 29cm	BR-FA0923	9F X 23cm	9F X 29cm
BR-RA0623	6F X 23cm	6F X 29cm	BR-FA1223	12F X 23cm	12F X 29cm

### 5. PRODUCT PERFORMANCE

Zero clearance is realized when using a combination of the Sheath catheter and dilator, which can effectively reduce the damage of blood vessel. The special design and excellent technology can effectively prevent blood backflow and air enter as well.

### 6. INSTRUCTION FOR USE

1. Remove the introducer sheath components from the packing box I and rinse the kit thoroughly with heparin saline.
2. Note: the intensity should be light and should be inserted from the center of the joint, otherwise it will damage the hemostatic valve that lead the sheath blood leak. Just as important, the dilator should be connected with the introducer sheath. Otherwise, the dilator will withdraw in the puncture of sheath kit and lead to the damage of the tip of the introducer sheath. It can cause damage to the blood vessels.
3. Insert the introducer sheath into the target vessel slowly along the guide wire.
4. The dilator and guide wire are slowly removed from the introducer sheath, but when the dilator is about to leave the introducer sheath, should quickly pull out the dilator, otherwise it will result in incomplete sealing and leakage of the hemostatic valve.

5. The catheter is inserted into the target vessel through the insertion of the sheath and is pushed into the desired location.
6. After the operation, the catheter is removed and then pulled out of the sheath and then compressed to stop bleeding after operation.
7. After use, destruction is performed according to scrap procedures for medical supplies.

7. **PRODUCT LABEL AND PACKAGING** (See the label or packaging)

1. Product name, model and specification.
2. Production date, batch lot.
3. Product license, registration number, product technology requirement.
4. Sterilization expiry date.

8. **CAUTION**

1. Procedure must be performed by trained medical personnel.
2. Sterilization is valid for three years, please confirm date of sterilization before use, do not use if expiry.
3. Do not use if luer taper leakage.
4. The device is designed and intended for **SINGLE USE ONLY**.
5. Do not resterilize and/or reuse, otherwise patients could get the risk of infection.
6. Sterilized by ethylene oxide. Do not use if the package damaged.



Caution



Use by date



Keep dry



Keep away from sunlight

REF

Catalogue number



Date of manufacture



Manufacturer



Do not use




Sterilized Using EO



Do not use if package is damaged



Batch code



BrosMed Medical B.V.

ADDRESS: Mgr. Buckxstraat 8, 6134 AP Sittard, the Netherlands

TEL: +31 850140901

Mail: [cs@brosmed.com](mailto:cs@brosmed.com)



Brun Medical Co., Ltd.

ADDRESS: Building B8, Entrepreneurship and Innovation Park, 88 Changtang Road,

Yunlong Demonstration Zone, Zhuzhou 412000 HUNAN, CHINA

TEL: +0731-22029666

WEB: [www.brunmd.com](http://www.brunmd.com)

Date of issue: 2022/05/05