

# Millau NC

0.035" system

## Non-Compliant OTW

Engineered for all PTA including AVF

Combination of strength and optimal control in routine and challenging cases


Designed to treat fibrotic and calcified lesions

Broad size range:

Diameter : 3.0 - 10.0mm

Length : 20 - 150mm

Five-fold(5 pleats)



# Specifications

Balloon Diameter (mm)	Balloon Length (mm)							Catheter Working Length (cm)	Sheath Compatibility (F)	Nominal Pressure (atm)	Rated Burst Pressure (atm)
	20	30	40	60	80	100	150				
3.0	V	V	V	V				40	5	12	22
4.0	V	V	V	V				40	5	12	22
5.0	V	V	V	V				40	5	12	20
6.0	V	V	V	V				40	6	12	20
7.0	V	V	V	V				40	6	12	18
8.0	V	V	V	V				40	6	12	18
3.0	V	V	V	V	V	V	V	70	5	12	22
4.0	V	V	V	V	V	V	V	70	5	12	22
5.0	V	V	V	V	V	V	V	70	5	12	20
6.0	V	V	V	V	V	V	V	70	6	12	20
7.0	V	V	V	V	V	V	V	70	6	12	18
8.0	V	V	V	V	V	V	V	70	6	12	18
9.0	V	V	V	V	V	V		70	7	12	18
10.0	V	V	V	V	V			70	7	12	18
3.0	V	V	V	V	V	V	V	90	5	12	22
4.0	V	V	V	V	V	V	V	90	5	12	22
5.0	V	V	V	V	V	V	V	90	5	12	20
6.0	V	V	V	V	V	V	V	90	6	12	20
7.0	V	V	V	V	V	V	V	90	6	12	18
8.0	V	V	V	V	V	V	V	90	6	12	18
9.0	V	V	V	V	V	V		90	7	12	18
10.0	V	V	V	V	V			90	7	12	18
3.0	V	V	V	V	V	V	V	150	5	12	22
4.0	V	V	V	V	V	V	V	150	5	12	22
5.0	V	V	V	V	V	V	V	150	5	12	20
6.0	V	V	V	V	V	V	V	150	6	12	20
7.0	V	V	V	V	V	V	V	150	6	12	18
8.0	V	V	V	V	V	V	V	150	6	12	18
9.0	V	V	V	V	V	V		150	7	12	18
10.0	V	V	V	V	V			150	7	12	18

Remark: Sizes within ■ have longer delivery time.



cNovate Medical BV  
 Address: Hogebrinkerweg 31, 3871 KM Hoevelaken, The Netherlands  
 Tel: +31 850140404 E-mail: cs@cNovate.eu Web: www.cNovate.eu

