

Skinny SC

0.018" system

Semi-Compliant OTW

Excellent for dilatation and lesion preparation

To treat a wide range of indications

Designed to treat resistant lesions

Consistent Rated Burst Pressure for all sizes

Tri-fold (3 pleats) for \varnothing 1.5-4.0mm

Five-fold(5 pleats) for \varnothing 4.5-10.0mm

Seamless tip-to-shaft design


Provides smooth transmission and maximized pushability

Low tip profile 0.020"

Facilitates lesion entry

hydrophilic coating

Latest hydrophilic coating technology to facilitate easy maneuvering through complex anatomy and lesions



Specifications

Balloon Diameter (mm)	Balloon Length (mm)												Catheter Working Length (cm)	Sheath Compatibility (F)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	
	5	10	15	20	30	40	60	80	100	120	150	200					
3.0				V	V	V	V							40	4	6	14
3.5				V	V	V	V							40	4	6	14
4.0				V	V	V	V							40	4	6	14
4.5				V	V	V	V							40	5	6	14
5.0				V	V	V	V							40	5	6	14
5.5				V	V	V	V							40	5	6	14
6.0				V	V	V	V							40	5	6	14
7.0				V	V	V	V							40	6	6	14
8.0				V	V	V	V							40	6	6	14
1.5	V	V	V	V	V	V	V	V	V	V	V	V	V	70	4	6	14
2.0	V	V	V	V	V	V	V	V	V	V	V	V	V	70	4	6	14
2.5	V	V	V	V	V	V	V	V	V	V	V	V	V	70	4	6	14
3.0	V	V	V	V	V	V	V	V	V	V	V	V	V	70	4	6	14
3.5	V	V	V	V	V	V	V	V	V	V	V	V	V	70	4	6	14
4.0	V	V	V	V	V	V	V	V	V	V	V	V	V	70	4	6	14
4.5	V	V	V	V	V	V	V	V	V	V	V	V	V	70	5	6	14
5.0	V	V	V	V	V	V	V	V	V	V	V	V	V	70	5	6	14
5.5	V	V	V	V	V	V	V	V	V	V	V	V	V	70	5	6	14
6.0	V	V	V	V	V	V	V	V	V	V	V	V	V	70	5	6	14
7.0	V	V	V	V	V	V	V	V	V	V	V	V	V	70	6	6	14
8.0	V	V	V	V	V	V	V	V	V	V	V	V		70	6	6	14
9.0	V	V	V	V	V	V	V	V	V	V	V			70	6	6	14
10.0	V	V	V	V	V	V	V	V						70	7	6	14
1.5	V	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
2.0	V	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
2.5	V	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
3.0	V	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
3.5	V	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
4.0	V	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
4.5	V	V	V	V	V	V	V	V	V	V	V	V	V	90	5	6	14
5.0	V	V	V	V	V	V	V	V	V	V	V	V	V	90	5	6	14
5.5	V	V	V	V	V	V	V	V	V	V	V	V	V	90	5	6	14
6.0	V	V	V	V	V	V	V	V	V	V	V	V	V	90	5	6	14
7.0	V	V	V	V	V	V	V	V	V	V	V	V	V	90	6	6	14
8.0	V	V	V	V	V	V	V	V	V	V	V	V		90	6	6	14
9.0	V	V	V	V	V	V	V	V	V	V	V			90	6	6	14
10.0	V	V	V	V	V	V	V	V						90	7	6	14
1.5	V	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
2.0	V	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
2.5	V	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
3.0	V	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
3.5	V	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
4.0	V	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
4.5	V	V	V	V	V	V	V	V	V	V	V	V	V	150	5	6	14
5.0	V	V	V	V	V	V	V	V	V	V	V	V	V	150	5	6	14
5.5	V	V	V	V	V	V	V	V	V	V	V	V	V	150	5	6	14
6.0	V	V	V	V	V	V	V	V	V	V	V	V	V	150	5	6	14
7.0	V	V	V	V	V	V	V	V	V	V	V	V	V	150	6	6	14
8.0	V	V	V	V	V	V	V	V	V	V	V	V		150	6	6	14
9.0	V	V	V	V	V	V	V	V	V	V	V			150	6	6	14
10.0	V	V	V	V	V	V	V	V						150	7	6	14

Remark: Sizes within ■ have longer delivery time.



cNovate Medical BV
 Address: Hogebrinkerweg 31, 3871 KM Hoevelaken, The Netherlands
 Tel: +31 850140404 E-mail: cs@cNovate.eu Web: www.cNovate.eu

