

# Zeeland SC

0.014" system

## Semi-Compliant OTW


Designed to treat BTK lesions



Excellent performance in routine  
and challenging cases

Designed to treat resistant lesions

For a wide range of indications



# Specifications

Balloon Diameter (mm)	Balloon Length (mm)												Catheter Working Length (cm)	Sheath Compatibility (F)	Nominal Pressure (atm)	Rated Burst Pressure (atm)
	5	10	15	20	30	40	60	80	100	120	150	200				
1.5	V	V	V	v	v	v	v	V	V	V	V	V	70	4	6	14
2.0	V	V	V	v	v	v	v	V	V	V	V	V	70	4	6	14
2.5	V	V	V	v	v	v	v	V	V	V	V	V	70	4	6	14
3.0	V	V	V	v	v	v	v	V	V	V	V	V	70	4	6	14
3.5	V	V	V	v	v	v	v	V	V	V	V	V	70	4	6	14
4.0	V	V	V	v	v	v	v	V	V	V	V	V	70	4	6	14
4.5	V	V	V	v	v	v	v	V	V	V	V	V	70	5	6	14
5.0	V	V	V	v	v	v	v	V	V	V	V	V	70	5	6	14
5.5	V	V	V	v	v	v	v	V	V	V	V	V	70	5	6	14
6.0	V	V	V	V	V	V	V	V	V	V	V	V	70	5	6	14
1.5	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
2.0	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
2.5	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
3.0	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
3.5	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
4.0	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
4.5	V	V	V	V	V	V	V	V	V	V	V	V	90	5	6	14
5.0	V	V	V	V	V	V	V	V	V	V	V	V	90	5	6	14
5.5	V	V	V	V	V	V	V	V	V	V	V	V	90	5	6	14
6.0	V	V	V	V	V	V	V	V	V	V	V	V	90	5	6	14
1.5	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
2.0	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
2.5	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
3.0	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
3.5	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
4.0	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
4.5	V	V	V	V	V	V	V	V	V	V	V	V	150	5	6	14
5.0	V	V	V	V	V	V	V	V	V	V	V	V	150	5	6	14
5.5	V	V	V	V	V	V	V	V	V	V	V	V	150	5	6	14
6.0	V	V	V	V	V	V	V	V	V	V	V	V	150	5	6	14

Remark: Sizes within  have longer delivery time.



cNovate Medical BV

Address: Hogebrinkweg 31, 3871 KM Hoevelaken, The Netherlands

Tel: +31 850140404 E-mail: cs@cNovate.eu Web: www.cNovate.eu



2797