

Millau NC

0.035" system

Non-Compliant OTW

Engineered for all PTA including AVF



Combination of strength and optimal control in routine and challenging cases


Designed to treat fibrotic and calcified lesions

Broad size range:

Diameter : 3.0 - 10.0mm Length : 20 - 150mm
Five-fold(5 pleats)

Excellent pushability

Seamless tip-to-shaft design for push force transmission



Lowest tip profile 0.038"

Facilitates lesion crossing with ease

Balanced non-compliant balloon

Controlled compliance to exert focal pressure on the lesion and safe dilatation

Specifications

Balloon Diameter (mm)	Balloon Length (mm)							Catheter Working Length (cm)	Sheath Compatibility (F)	Nominal Pressure (atm)	Rated Burst Pressure (atm)
	20	30	40	60	80	100	150				
3.0	V	V	V	V				40	5	12	22
4.0	V	V	V	V				40	5	12	22
5.0	V	V	V	V				40	5	12	20
6.0	V	V	V	V				40	6	12	20
7.0	V	V	V	V				40	6	12	18
8.0	V	V	V	V				40	6	12	18
3.0	V	V	V	V	V	V	V	70	5	12	22
4.0	V	V	V	V	V	V	V	70	5	12	22
5.0	V	V	V	V	V	V	V	70	5	12	20
6.0	V	V	V	V	V	V	V	70	6	12	20
7.0	V	V	V	V	V	V	V	70	6	12	18
8.0	V	V	V	V	V	V	V	70	6	12	18
9.0	V	V	V	V	V	V		70	7	12	18
10.0	V	V	V	V	V			70	7	12	18
3.0	V	V	V	V	V	V	V	90	5	12	22
4.0	V	V	V	V	V	V	V	90	5	12	22
5.0	V	V	V	V	V	V	V	90	5	12	20
6.0	V	V	V	V	V	V	V	90	6	12	20
7.0	V	V	V	V	V	V	V	90	6	12	18
8.0	V	V	V	V	V	V	V	90	6	12	18
9.0	V	V	V	V	V	V		90	7	12	18
10.0	V	V	V	V	V			90	7	12	18
3.0	V	V	V	V	V	V	V	150	5	12	22
4.0	V	V	V	V	V	V	V	150	5	12	22
5.0	V	V	V	V	V	V	V	150	5	12	20
6.0	V	V	V	V	V	V	V	150	6	12	20
7.0	V	V	V	V	V	V	V	150	6	12	18
8.0	V	V	V	V	V	V	V	150	6	12	18
9.0	V	V	V	V	V	V		150	7	12	18
10.0	V	V	V	V	V			150	7	12	18

Remark: Sizes within ■ have longer delivery time.



cNovate Medical BV

Address: Terminalweg 15, 3821 AJ Amersfoort, The Netherlands

Tel: +31 850140404 E-mail: cs@cNovate.eu Web:www.cNovate.eu

