



Trift sc

0.035" system

Semi-Compliant OTW
Engineered for all PTA including AVF

Combines strength and optimal control in routine and challenging cases


To treat fibrotic and calcified lesions

Broad size range:

Diameter : 3.0 - 10.0mm

Length : 20 - 200mm

Five-fold(5 pleats)



Excellent pushability

Seamless tip-to-shaft design for smooth push force transmission

Lowest tip profile 0.038"

Facilitates lesion crossing with ease

hydrophilic coating

Latest hydrophilic coating technology to facility easy maneuvering through complex anatomy and lesions

Specifications

Balloon Diameter (mm)	Balloon Length (mm)									Catheter Working Length (cm)	Sheath Compatibility (F)	Nominal Pressure (atm)	Rated Burst Pressure (atm)
	20	30	40	60	80	100	120	150	200				
3.0	V	V	V	V						40	5	6	14
4.0	V	V	V	V						40	5	6	14
5.0	V	V	V	V						40	5	6	14
6.0	V	V	V	V						40	6	6	14
7.0	V	V	V	V						40	6	6	14
8.0	V	V	V	V						40	6	6	14
3.0	V	V	V	V	V	V	V	V	V	70	5	6	14
4.0	V	V	V	V	V	V	V	V	V	70	5	6	14
5.0	V	V	V	V	V	V	V	V	V	70	5	6	14
6.0	V	V	V	V	V	V	V	V	V	70	6	6	14
7.0	V	V	V	V	V	V	V	V	V	70	6	6	14
8.0	V	V	V	V	V	V	V	V		70	6	6	14
9.0	V	V	V	V	V	V	V			70	7	6	14
10.0	V	V	V	V	V					70	7	6	14
3.0	V	V	V	V	V	V	V	V	V	90	5	6	14
4.0	V	V	V	V	V	V	V	V	V	90	5	6	14
5.0	V	V	V	V	V	V	V	V	V	90	5	6	14
6.0	V	V	V	V	V	V	V	V	V	90	6	6	14
7.0	V	V	V	V	V	V	V	V	V	90	6	6	14
8.0	V	V	V	V	V	V	V	V		90	6	6	14
9.0	V	V	V	V	V	V	V			90	7	6	14
10.0	V	V	V	V	V					90	7	6	14
3.0	V	V	V	V	V	V	V	V	V	150	5	6	14
4.0	V	V	V	V	V	V	V	V	V	150	5	6	14
5.0	V	V	V	V	V	V	V	V	V	150	5	6	14
6.0	V	V	V	V	V	V	V	V	V	150	6	6	14
7.0	V	V	V	V	V	V	V	V	V	150	6	6	14
8.0	V	V	V	V	V	V	V	V		150	6	6	14
9.0	V	V	V	V	V	V	V			150	7	6	14
10.0	V	V	V	V	V					150	7	6	14

Remark: Sizes within ■ have longer delivery time.



cNovate Medical BV
 Address: Hogebrinkerweg 31, 3871 KM Hoevelaken, The Netherlands
 Tel: +31 850140404 E-mail: cs@cNovate.eu Web: www.cNovate.eu

