

Trift sc

0.035" system

Semi-Compliant OTW
Engineered for all PTA including AVF

Combines strength and optimal control in routine and challenging cases

To treat fibrotic and calcified lesions

Broad size range:

Diameter : 3.0 - 10.0mm

Length : 20 - 200mm

Five-fold(5 pleats)

Excellent pushability


Seamless tip-to-shaft design for smooth push force transmission

Lowest tip profile 0.038"

Facilitates lesion crossing with ease

hydrophilic coating

Latest hydrophilic coating technology to facility easy maneuvering through complex anatomy and lesions



Specifications

Balloon Diameter (mm)	Balloon Length (mm)									Catheter Working Length (cm)	Sheath Compatibility (F)	Nominal Pressure (atm)	Rated Burst Pressure (atm)
	20	30	40	60	80	100	120	150	200				
3.0	V	V	V	V						40	5	6	14
4.0	V	V	V	V						40	5	6	14
5.0	V	V	V	V						40	5	6	14
6.0	V	V	V	V						40	6	6	14
7.0	V	V	V	V						40	6	6	14
8.0	V	V	V	V						40	6	6	14
3.0	V	V	V	V	V	V	V	V	V	70	5	6	14
4.0	V	V	V	V	V	V	V	V	V	70	5	6	14
5.0	V	V	V	V	V	V	V	V	V	70	5	6	14
6.0	V	V	V	V	V	V	V	V	V	70	6	6	14
7.0	V	V	V	V	V	V	V	V	V	70	6	6	14
8.0	V	V	V	V	V	V	V	V		70	6	6	14
9.0	V	V	V	V	V	V	V			70	7	6	14
10.0	V	V	V	V	V					70	7	6	14
3.0	V	V	V	V	V	V	V	V	V	90	5	6	14
4.0	V	V	V	V	V	V	V	V	V	90	5	6	14
5.0	V	V	V	V	V	V	V	V	V	90	5	6	14
6.0	V	V	V	V	V	V	V	V	V	90	6	6	14
7.0	V	V	V	V	V	V	V	V	V	90	6	6	14
8.0	V	V	V	V	V	V	V	V		90	6	6	14
9.0	V	V	V	V	V	V	V			90	7	6	14
10.0	V	V	V	V	V					90	7	6	14
3.0	V	V	V	V	V	V	V	V	V	150	5	6	14
4.0	V	V	V	V	V	V	V	V	V	150	5	6	14
5.0	V	V	V	V	V	V	V	V	V	150	5	6	14
6.0	V	V	V	V	V	V	V	V	V	150	6	6	14
7.0	V	V	V	V	V	V	V	V	V	150	6	6	14
8.0	V	V	V	V	V	V	V	V		150	6	6	14
9.0	V	V	V	V	V	V	V			150	7	6	14
10.0	V	V	V	V	V					150	7	6	14

Remark: Sizes within have longer delivery time.



cNovate Medical BV

Address: Terminalweg 15, 3821AJ Amersfoort, The Netherlands

Tel: +31 850140404 E-mail: cs@cNovate.eu Web: www.cNovate.eu



2797